EN

Important information - please read prior to use

Guided Growth System (Bone Plate)

Orthofix Srl
Via delle Nazioni 9 - 37012 Bussolengo (VR) Italy
Telephone +39 045 671 9000 - Fax +39 045 671 9380

INDICATIONS FOR USE

The Guided Growth System is designed for the express and sole purpose of redirecting the angle of growth of long bone(s). This is useful for gradually correcting angular growth deformities in growing children. Specific conditions/diseases for which the device will be indicated include: valgus, varus or flexion, extension deformities of the knee (femur and/or tibia); valgus, varus, or planar flexion deformities of the ankle; valgus or varus deformities of the elbow (humerus); and radial or ulnar deviation, flexion or extension deformities of the wrist (radius).

DESCRIPTION

The Guided Growth System is designed for the gradual correction of pediatric deformities in both the upper and lower extremities. The device can be used for correction of congenital deformities as well as correction of acquired deformities, provided that the physis (growth plates) are not fused.

The plate features a contoured waist and low profile for pediatric usage. There is a center hole in the plate for a temporary guide pin to be implanted to ensure accurate application of the plate. The Plate is attached to the external surface of the bone over the growth plate by screws. These screws are not locked to the plate, but rather are allowed to swivel and diverge in their position as bone growth occurs. The implant acts like a flexible hinge, permitting growth at the growth plate to gradually straighten the limb. Immediately after implantation, the patient is allowed mobility and weight bearing.

CONTRAINDICATIONS

- The Guided Growth System should not be used for adult deformities; skeletal maturity precludes further growth.
- The Guided Growth System should not be used where the growth plate has closed due to trauma or infection since further growth is precluded.

WARNINGS

- The stainless steel plates may be bent to a 20 degree offset for more severe anatomical bone contours. Bending beyond 20 degree offset could result in injury or reoperation due to breakage.
- Particular care should be taken that bone screws do not enter the joints or damage the growth plates in children.
- Bone plates and screws must not be reused. If any implants has come into contact with any body fluid it should be considered to have been used.
- The Guided Growth System has not been evaluated for safety and compatibility in the MRI environment. The Guided Growth System has not been tested for heating or migration in the MR environment.
- This device is not approved for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic or lumbar spine.

PRECAUTIONS

- It is essential that proper operative technique be followed for implantation. Advance bone screws in an alternating manner until the screws are fully seated in the plate and the plate is flush with the bone. Ensure positioning after insertion with fluoroscopy. Please refer to the appropriate Operative Technique.
- Examine all components carefully prior to use. Product integrity, sterility (in the case of sterile products) and performance are assured only if the packaging is undamaged. DO NOT USE if packaging is compromised or if a component is believed to be faulty, damaged or suspect.
- Do not combine Guided Growth System implantable components with those from other manufacturers.
- Combining implants of dissimilar metals may result in galvanic corrosion.

POSSIBLE ADVERSE EFFECTS

- Damage to the growth plate can occur if the screws enter the growth plate. It is important to follow the recommended implantation technique to avoid this.
- Failure to achieve the desired angular correction or possible over-correction can occur.
- Breaking of the plate or loosening, bending, or breaking of the bone screws.
- Re-operation may be necessary to replace or remove the bone plate and screws.
- Complications associated with metal sensitivity.
- Intrinsic risks associated with anesthesia and surgery.

IMPORTANT

A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons or device failure that require surgical re-intervention to remove or replace the plate or bone screws. Proper implantation and patient selection will greatly affect the results. Patients must be closely followed since achievement of the correct angle may require plate removal as early as 3 to 12 months after insertion. Rebound overgrowth after plate removal can sometimes require repeat surgery.

MATERIALS

The Guided Growth Plates and Bone Screws are manufactured in 2 different configurations:

- made from implant grade 316L stainless steel;
 - made from implant grade titanium alloy.
- Labelling and product codes permit an easy identification of configurations.

STERILITY

Orthofix provides Guided Growth Plates and Screws in sterile or non-sterile configuration. Labelling and product codes permit an easy identification of configurations. The products supplied NON-STERILE require cleaning and sterilization prior to use, following the respective recommended procedures.

Cleaning, Sterilization and Maintenance

Cleaning is an essential pre-requisite to ensure effective disinfection or sterilization. All instruments must be cleaned prior to use and after each re-use of the instruments. Cleaning may be done following the validated cleaning process described below and using washing instruments and machines validated by the hospital.

WARNINGS

- Aluminum based instruments are damaged by alkaline (pH>7) detergents and solutions.
- Anodised coating is damaged by detergents with free halogen ions or sodium hydroxide.
- Detergents and disinfectants with fluoride, chloride, bromide, iodide or hydroxy ions MUST NOT be used.
- PRODUCTS LABELED FOR SINGLE-USE MUST NOT BE REUSED.

Limitations and restrictions on reprocessing

- Repeated reprocessing has minimal effect on reusable instruments;
- End of life is normally determined by wear and damage due to use.

Point of use

- It is recommended that instruments are reprocessed as soon as is reasonably practicable following use, putting them inside a sink filled with cold water (< 40°C) for at least 10 minutes and removing gross soiling with a soft cloth or a soft brush.
- Don't use a floating detergent or hot water as this can cause the fixation of residue which may influence the result of the reprocessing process.

Containment and transportation

- Follow hospital protocols when handling contaminated and bio-hazardous materials;
- All used surgical instruments should be regarded as contaminated. Their handling, collection and transportation should be strictly controlled to minimize any possible risks to patients, personnel and any area of the healthcare facility.

- Preparation for decontamination**
- Disassemble the devices where applicable. See Orthofix operative techniques for further details.

Cleaning: manual

- Soak the single components in the cleaning solution. Orthofix recommends the use of a pH-neutral enzymatic cleaning agent. Please refer to detergent's manufacturer datasheet for concentration of the solution, required time and temperature.
- Brush the single components in the cleaning solution with a surface soft brush until all visible soiling is removed. Use a soft brush to remove residuals from lumens
- Ultrasonic clean the single components in a degassed cleaning solution. Orthofix recommends the use of a pH-neutral or slightly alkaline solution. Please refer to detergent's manufacturer datasheet for concentration of the solution, required time and temperature.
- Brush the single components in running tap water
- Carefully hand-dry using absorbent, non-shedding cloth or industrial dryer.
- Rinse the components in sterile or deionized water
- Carefully hand-dry using absorbent, non-shedding cloth or industrial dryer.

Cleaning: automated

When the devices to be cleaned have lumens or present complexity, a preliminary manual cleaning may be required.

- Place all the instruments into washer baskets.
- Orient instruments into the automated washer's carriers as recommended by the washer manufacturer.
- Orthofix recommends the use of a pH-neutral or slightly alkaline solution. Please refer to detergent's manufacturer datasheet for concentration of the solution, required time and temperature.

Orthofix recommends that cycle steps are at least as follows:

- prewash at low temperatures
- main wash
- rinse with demineralized water
- thermal disinfection
- drying
- Rinse the components in sterile or deionized water
- Carefully hand-dry using absorbent, non-shedding cloth or industrial dryer.

Maintenance, inspection and testing

All instruments and product components should be visually inspected for cleanliness.

If the equipment does not appear visually clean, cleaning steps must be repeated until the products do not result visually clean.

All instruments and product components should be visually inspected for any signs of deterioration that may cause failure in use (such as cracks or damage to surfaces) and functions tested before being sterilized(see detailed Operative Technique Manuals and Instructions for use). Particular attention should be given to:

- Cannulated devices (NB: cannulated drill-bits are single-patient use only).
- Cutting edges. Discard blunt or damaged instruments.

- If a component or instrument is believed to be faulty, damaged or suspect, it should NOT BE USED.

- When instruments form part of an assembly, check assembly with matching components.

Packaging

- Wrap the tray before sterilization with an approved sterilization wrap or insert it into a rigid sterilization container in order to prevent contamination after sterilization.
- Do not re-include additional systems or instruments in the sterilization tray. Sterility cannot be guaranteed if the sterilization tray is overloaded.

Sterilization

- Steam sterilization is recommended. Gas plasma and EtO sterilization should be avoided as they were not validated for Orthofix products;
- Use a validated, properly maintained and calibrated steam sterilizer;
- Do not exceed 140°C (284°F);
- Do not stack trays during sterilization
- Sterilize by steam autoclaving, utilizing a fractioned pre-vacuum cycle or gravity cycle according to the table below:

STEAM STERILIZER TYPE	GRAVITY	PRE-VACUUM
Preconditioning pulses	N/A	4 pulses
Exposure Temperature	132°C (270°F)	132°C (270°F)
Minimum Holding Time	15 minutes	4 minutes
Drying Time	30 minutes	30 minutes

Storage

Store the sterilized instrument in a dry, clean, and dust-free environment at room temperature.

SERVICE

For any information please contact your local distributor.

CAUTION

Federal law (U.S.) restricts this device to sale by or on the order of a physician.

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