THE GOTFRIED PERCUTANEOUS COMPRESSION PLATE (PC.C.P.) 3	EN	GOTFRIED PERKUTAN KOMPRESSIONSPLATTA (PC.C.P.) 35	SV
PLACCA DI COMPRESSIONE PERCUTANEA DI GOTFRIED (PC.C.P.) 6	IT	Η Διαδερμική Πλάκα Συμπίεσης Gotfried (PC.C.P.) 38	EL
PLAQUE DE COMPRESSION PERCUTANÉE (PC.P.C.) DE GOTFRIED 9	FR	ゴットフリード経皮的 コンプレッションプレート(PC.C.P.) 42	JP
DIE PERKUTANE KOMPRESSIONSPLATTE NACH GOTFRIED (PC.C.P.) 13	DE	GOTFRIED 经皮加压钢板 (PC.C.P.) 45	CN
PLACA DE COMPRESIÓN PERCUTÁNEA GOTFRIED (PC.C.P.) 16	ES	GOTFRIEDOVA PERKUTÁNNÍ KOMPRESNÍ DLAHA (PC.C.P.) 48	CS
GOTFRIED PERKUTAN KOMPRESSIONSSKINNE (PC.C.P.) 19	DA	PLYTKA GOTFRIEDA DO KOMPRESJI PRZEZSKÓRNEJ (PC.C.P.) 51	PL
GOTFRIED PERKUTAANINEN KOMPRESSIOLEVY (PC.C.P.) 22	FI	GOTFRIED PERKUTANA KOMPRESIJSKA PLOŠČA (PC.C.P.) 54	SL
PERKUTAN GOTFRIED KOMPRESJONSPLATE (PC.C.P.) 25	NO	GOTFRIED (PC.C.P.) 58	KR
DE GOTFRIED PERCUTANE COMPRESSIEPLAAT (PC.C.P.) 28	NL	 ORTHOFIX® ORTHOFIX Srl Via Delle Nazioni 9 - 37012 Bussolengo (Verona) - Italy Telephone 0039-0456719000 - Fax 0039-0456719380	
PLACA DE COMPRESSÃO PERCUTÁNEA GOTFRIED (PC.C.P.) 32	PT		

THE GOTFRIED PERCUTANEOUS COMPRESSION PLATE (PC.C.P.)

DESCRIPTION & INDICATIONS FOR USE

The Gotfried Percutaneous Compression Plate (PC.C.P.) for trochanteric femoral fractures is intended as a means of bone stabilization in the management of fractures. The use of external supports (e.g. walking aids) is recommended as a part of the treatment.

The Gotfried PC.C.P. is designed specifically for the fixation of acute trochanteric fractures of the proximal femur. It is not suitable for the fixation of intra-capsular fractures of the femoral neck, or for fractures extending into the subtrochanteric region.

The product is intended for professional use only. Surgeons who supervise the use of the product must have full awareness of orthopaedic fixation procedures as well as adequate understanding of the philosophy of the product. To promote the proper use of the product, and establish an effective promotional and training tool, Orthofix has developed an operative technique manual as well as an operative video, which is available as a complimentary service for surgeons who have adopted the Gotfried PC.C.P.

If you wish to receive a personal copy of this manual, please contact Orthofix or its local authorized representative.

CONTRAINDICATIONS

The Gotfried PC.C.P. is not designed or sold for any use except as indicated.

Use of the Gotfried PC.C.P. is contraindicated in the following situations:

1. Active infection.
2. Acute medical conditions including: impaired blood supply, pulmonary insufficiency (i.e. ARDS, fat embolism) or latent infection.
3. Foreign body sensitivity. Where material sensitivity is suspected, tests should be made prior to implantation.

WARNINGS & PRECAUTIONS

1. The Gotfried Percutaneous Compression Plate (PC.C.P.) is for the treatment of acute trochanteric femoral fractures only. It must not be used for intracapsular femoral neck fractures or for subtrochanteric fractures. It should only be used if there is a proper fracture table available, so that the fracture can be held in a stable reduced position. The use of a Posterior Reduction Device (PORD) is required to achieve this. It is very important that the operating surgeon has received instruction on the insertion technique for the PC.C.P., and that the full set of instruments is available. The full range of femoral neck screws must be available so that the correct size can be selected.
2. The Gotfried PC.C.P. should not be used if the fracture is open, if the overlying skin is damaged, or if active infection is present. In this case External Fixation is preferable.
3. For a successful result after insertion of a PC.C.P. the patient should be mobilised immediately, weightbearing as much as can

- be tolerated. Full weightbearing is encouraged from the beginning.
4. Fracture distraction for any time should be avoided during the operation, and fractures should never be locked in distraction, since this may cause delayed union.
 5. Correct handling of the instrumentation and implant is extremely important. Surgical theatre staff should avoid any notching or scratching of the implants, as these factors may produce internal stresses which may become a focus for possible implant breakage.
 6. Any fixation device may break if subjected to the increased loading caused by delayed or non-union.
 7. Implants and locking screws must never be reused.
 8. Patients should be instructed to report any adverse or unanticipated effects to the treating surgeon.
 9. Nuclear Magnetic Resonance imaging should not be used in any segment to which an implant is applied.
 10. All equipment should be carefully examined prior to use to assure proper working condition. If a component or instrument is believed to be faulty, damaged or suspect, it should NOT BE USED.
 11. Additional equipment may be required for application and removal such as soft tissue retractors, etc.

POSSIBLE ADVERSE EFFECTS

1. Non-union or delayed union, which may lead to implant breakage.
2. Metal sensitivity, or allergic reaction to a foreign body.
3. Pain, discomfort, or abnormal sensations due to the presence of the device.

Warning: This device is not approved for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic or lumbar spine.

IMPORTANT

A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons or device failure which require surgical reintervention to remove or replace the internal fixation device.

Preoperative and operative procedures including knowledge of surgical techniques and proper selection and placement of the device are important considerations in the successful utilization of the devices by the surgeon.

Proper patient selection and the patient's ability to comply with physician instructions and follow prescribed treatment regimen will greatly affect the results. It is important to screen patients and select optimal therapy given physical and/or mental activity requirements and/or limitations. If a surgical candidate exhibits any contraindication or is predisposed to any contraindication, DO NOT USE the PC.C.P.

MATERIALS

The Gotfried PC.C.P. System is comprised of stainless steel and aluminium alloy components. Those components which come into contact with patients are made of surgical grade stainless steel.

STERILE & NON-STERILE PRODUCT

Orthofix provides certain devices STERILE while others are provided NON-STERILE.

Please review the product label to determine the sterility of each device.

STERILE

Devices or kits provided STERILE are labeled as such.

Contents of package are STERILE unless package is opened or damaged. Do not use if package is opened or damaged.

NON-STERILE

All components are provided NON-STERILE. The manufacturer recommends that all NON-STERILE components be properly cleaned and sterilized following the recommended cleaning and sterilization procedures.

Product integrity and performance are assured only if the packaging is undamaged.

CLEANING & MAINTENANCE (Attention, see instructions PQ ISP)

Prior to use, NON-STERILE product must be cleaned using a mixture of 70% medical grade alcohol and 30% distilled water. After cleaning, the device and/or system components should be thoroughly rinsed in sterile distilled water and dried using clean non-woven fabric. Lubricate all parts with Orthofix silicone oil, which is available through Orthofix or your local representative. If the product is of the reusable type, after it has been removed from the patient, it should be disassembled completely, all components cleaned using either 12 vol. hydrogen peroxide or a recommended detergent. To prevent corrosion, the components must be kept dry, and detergents with fluoride, chloride, bromide, iodide or hydroxyl ions must be avoided when cleaning, as they will damage the black anodized coating on any Orthofix products, and this may initiate the process of stress corrosion. Before sterilization, all components should be inspected, since damage to the surface of metal components can reduce strength and fatigue resistance, and may lead to corrosion. If components are damaged in any way, they should be exchanged immediately for new ones. Assembly of the product/system should then be carried out to ensure that all components are present.

ANY DEVICE WHICH IS LABELLED "SINGLE USE ONLY" MUST NEVER BE REUSED. ORTHOFIX IS ONLY RESPONSIBLE FOR SAFETY AND EFFECTIVENESS FOR THE FIRST PATIENT USE OF SINGLE USE DEVICES. The institution or practitioner bears full responsibility for any subsequent use of these devices.

STERILIZATION

The recommended, validated sterilization cycle is:

Method	Cycle	Temperature	Exposure Time
Steam	Pre-Vacuum	132° - 135°C [270° - 275°F]	Minimum 10 minutes

Sterility cannot be assured if the sterilization tray is overloaded. Do not overload the sterilization tray or include additional implants or instruments from any source. Orthofix is only responsible for safety and effectiveness for the initial patient use of the fixation devices. The institution or practitioner bears full responsibility for any subsequent use of the devices.

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.

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. ORTHOFIX

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: ()

All Orthofix internal and external fixation products should be used together with their corresponding Orthofix implants, components and accessories. Their application should be performed with the specific Orthofix instrumentation, following carefully the surgical technique recommended by the manufacturer in the appropriate Operative Technique Manual.

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