ORTHOFIX® PREFIX ² 2	EN	ORTHOFIX® PREFIX ² 42	SV
ORTHOFIX® PREFIX ² 6	IT	ORTHOFIX® PREFIX ² 46	EL
ORTHOFIX® PREFIX ² 10	FR	ORTHOFIX® PREFIX ² 50	JP
ORTHOFIX® PREFIX ² 14	DE	ORTHOFIX® PREFIX ² 54	CN
ORTHOFIX® PREFIX ² 18	ES	ORTHOFIX® PREFIX ² 58	CS
ORTHOFIX® PREFIX ² 22	DA	ORTHOFIX® PREFIX ² 62	PL
ORTHOFIX® PREFIX ² 26	FI	ORTHOFIX® PREFIX ² 66	SL
ORTHOFIX® PREFIX ² 30	NO	ORTHOFIX® PREFIX ² 70	KR
ORTHOFIX® PREFIX ² 34	NL		
ORTHOFIX® PREFIX ² 38	PT		

ORTHOFIX Srl
Via Delle Nazioni 9 - 37012 Bussolengo (Verona) - Italy
 Telephone 0039-0456719000 - Fax 0039-0456719380

ORTHOFIX® PREFIX²

GENERAL WARNINGS

- In order to counter-act the tightening forces it is advisable to use the Clamp Holder (92017) to grip the Rod Clip Clamp, the Screw Clip Clamp or the Transfixing Pin Clamp while locking them with the Universal T-Wrench or the 5 mm Allen Wrench.
- When using the Multi-Screw Clamp, two rods must always be used to increase frame stiffness.

ASPECTS OF MRI COMPATIBILITY

MRI machine: PHILIPS Achieva 1.5T A series

Device Position: isocenter of MRI scanner

Magnetic Field Strength: 1.5 Tesla

Spatial gradient field: 9.37mT/m

Please note that temperature changes reported apply to the designated MRI system and characteristics used. If a different MRI system is used, temperature changes may vary. Contact is advised with the supplier of the MRI system and Orthofix.

Prefix² System is designed as MRI conditional up to 1.5 Tesla. The MRI components have been tested according to ASTM Standards F2052, F2182, F2213 and F2119. All Prefix² Clamps are blue and MRI conditional. They are marked with "MR". Prefix² System can only be guaranteed for MRI when using Orthofix XCaliber Bone Screws and the following rods to build a frame.

CLAMPS

Code	Description
92012	Blue Rod Clip Clamp
92016	Blue Screw Clip Clamp
92026	Blue Transfixing Clip Clamp
92031	Blue Multi Screw Clamp

RODS

Code	Description
92125	Rod 125 mm long
92175	Rod 175 mm long
92275	Rod 275 mm long
92350	Rod 350 mm long

XCALIBER SCREWS

Total Length (mm)	Thread Length (mm)						
	30	40	50	60	70	80	90
150	911530	911540	911550	911560	911570	911580	911590
260	912630	912640	912650	912660	912670	912680	912690

It has been shown by specific MRI tests that the Prefix[®] System may be used for patients undergoing MRI procedures using up to 1.5 Tesla MR systems if certain specific conditions are followed. Four commonly used frames have been tested for MRI use at 1.5 Tesla. The results are as follows:

Blue Screw Clip Clamp in Single Bar Configuration

- A maximum temperature increase of 6.0 °C has been measured at the tip of the most external screw (Maximum Whole Body averaged SAR of 2W/kg for 3.5 minutes of scanning)

Code	Description
1x92350	Rod 350 mm long
4x911560	Self-Drilling XCaliber Screw, L. 150 mm, thread length 60 mm
4x92016	Blue Screw Clip Clamp

Blue Multi Screw Clamp

- A maximum temperature increase of 4.0 °C has been measured at the tip of the indicated screws (Maximum Whole Body averaged SAR of 2W/kg for 6 minutes of scanning)

Code	Description
2x92350	Rod 350 mm long
4x911560	Self-Drilling XCaliber Screw, L. 150 mm, thread length 60 mm
4x92012	Blue Rod Clip Clamp
2x92031	Blue Multi Screw Clamp

Z Configuration


- A maximum temperature increase of 7.0 °C has been measured at the tip of the most external screw (Maximum Whole Body averaged SAR of 2W/kg for 6 minutes of scanning)

Code	Description
2x92125	Rod 125 mm long
1x92275	Rod 275 mm long
2x911560	Self-Drilling XCaliber Screw, L. 150 mm, thread length 60 mm
2x911540	Self-Drilling XCaliber Screw, L. 150 mm, thread length 40 mm
4x92016	Blue Screw Clip Clamp
2x92012	Blue Rod Clip Clamp

Ankle Configuration

- A maximum temperature increase of 2.1 °C has been measured at the tip of the indicated screw (Maximum Whole Body averaged SAR of 2W/kg for 6 minutes of scanning)
- At the indicated clamp a temperature increase of 12.2 °C (with reference to room temperature) has been measured

Code	Description
1x92125	Rod 125 mm long
1x92175	Rod 175 mm long
1x92275	Rod 275 mm long
2x911540	Self-Drilling XCaliber Screw, L. 150 mm, thread length 40 mm
1x92080	Transfixing Pin, thread length 80 mm, thread Ø 5 mm, shaft Ø 4 mm
2x92016	Blue Screw Clip Clamp
2x92012	Blue Rod Clip Clamp
2x92026	Blue Transfixing Pin Clamp

These tests have been performed in areas where the greatest temperature increase is expected with commonly used frames: at the tip and surrounding the tip of the most external screw. Due to the versatility of the system, an unlimited number of frames can be built which makes it impossible to test every construct.

There are factors that can influence these results (e.g. the number of screws used). Therefore, it is recommended that each frame be evaluated by a radiologist or MR scientist before the MRI procedure is undertaken to ensure patient safety. Since different frame configurations and frame sizes might lead to higher temperature increases, Orthofix recommends that the SAR settings are minimised as much as possible.

None of the components should move or migrate in the 1.5 Tesla MRI environments. Non-clinical testing has not been performed to rule out the possibility of component movement or migration at static magnetic field strengths higher than 1.5 Tesla or maximum spatial gradients higher than 9.37mT/m.

IMAGE ARTIFACTS

MR image quality may be compromised if the area of interest is in the same vicinity or relatively close to the position of the device. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this implant. Representative components used to assemble Pre-fix® have been evaluated in MRI chamber. Typical values of artifacts sizes are reported below.

Code	Description	
1x92275	Rod 275 mm long	
1x92012	Blue Rod Clip Clamp	
1x92016	Blue Screw Clip Clamp	
1x911560	Self-Drilling XCaliber Screw, L. 150 mm, thread length 60 mm	

Pulse Sequence	Spin-echo	Gradient-echo
TR (ms)	500	334
TE (ms)	4,8	4,6
Slice Thickness (mm)	5	5
Matrix size	230x288	230x288
Flip Angle	90°	30°
Imaging Plane	perpendicular to static field	perpendicular to static field
Artifact Size (mm)	16	34

Code	Description		
1x92350	Rod 350 mm long		
1x92012	Blue Rod Clip Clamp		
1x92031	Blue Multi Screw Clamp		
2x911560	Self-Drilling XCaliber Screw, L. 150 mm, thread length 60 mm		
Pulse Sequence	Spin-echo	Gradient-echo	
TR (ms)	500	334	
TE (ms)	4,8	4,6	
Slice Thickness (mm)	5	5	
Matrix size	230x288	230x289	
Flip Angle	90°	30°	
Imaging Plane	perpendicular to static field	perpendicular to static field	
Artifact Size (mm)	20	46	

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All Orthofix internal and external fixation products should be used together with their corresponding Orthofix implants, components and accessories. Their application should be performed with the specific Orthofix instrumentation, following carefully the surgical technique recommended by the manufacturer in the appropriate Operative Technique Manual.

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